

Statement of Deficiencies	(X1) Provider/Supplier/CLIA Identification Number 29D2017786	(X3) Date Survey Completed 12/11/2018
Name of Provider or Supplier Las Vegas Skin & Cancer Clinics-Centennial Hills	Street Address, City, State 9780 W Skye Canyon Park Dr Ste 100, Las Vegas, NV	
For information on the provider's plan to correct this deficiency, please contact the provider or the state survey agency.		

(X4) ID Prefix Tag	Summary Statement of Deficiencies
D0000	<p>This Statement of Deficiencies was created as a result of an on-site CLIA recertification survey conducted at your facility on December 11, 2018. The findings and conclusions of any investigation by the Division of Public and Behavioral Health shall not be construed as prohibiting any criminal or civil investigations, actions or other claims for relief that may be available to any party under applicable federal, state, or local laws.</p>
D5785	<p>CORRECTIVE ACTIONS CFR(s): 493.1282(b)(3)</p> <p>(b) The laboratory must document all corrective actions taken, including actions taken when any of the following occur: (b)(3) The criteria for proper storage of reagents and specimens, as specified under 493.1252(b), are not met.</p> <p>This STANDARD is not met as evidenced by: Based on a review of temperature recording logs from January 1, 2017 through December 11, 2018 and an interview with the director of compliance, the laboratory failed to document daily temperatures to monitor the proper storage of reagents and specimens. Findings include: 1. The laboratory failed to document daily temperatures for the cryostat machine of minus 20 to minus 30 degrees Centigrade, the laboratory refrigerator of 2 to 8 degrees Centigrade and the laboratory room temperature of 22 to 25 degrees Centigrade. 2. There were no temperature recordings for 21 of 23 days for the cryostat machine, the laboratory refrigerator and the laboratory room temperature for the month of March 2017. 3. There were no temperature recordings for 9 of 23 days for the cryostat machine, the laboratory refrigerator and the laboratory room temperature for the month of June 2017. 4. There were no temperature recordings for 15 of 22 days for the cryostat machine, the laboratory refrigerator and the laboratory</p>

room temperature for the month of December 2017. This was confirmed by the director of compliance on December 11, 2018 at approximately 10:30 am. The laboratory performs approximately 250 histopathology patient tests annually.