

Statement of Deficiencies	(X1) Provider/Supplier/CLIA Identification Number 44D2147090	(X3) Date Survey Completed 08/19/2026
Name of Provider or Supplier Rhsc Medical Center Of Jonesborough	Street Address, City, State 600 Boones Creek Rd, Jonesborough, TN	
For information on the provider's plan to correct this deficiency, please contact the provider or the state survey agency.		

(X4) ID Prefix Tag	Summary Statement of Deficiencies
D5415	TEST SYSTEMS, EQUIPMENT, INSTRUMENTS, REAGENT CFR(s): 493.1252(c) (c) Reagents, solutions, culture media, control materials, calibration materials, and other supplies, as appropriate, must be labeled to indicate the following: (c)(1) Identity and when significant, titer, strength or concentration. (c)(2) Storage requirements. (c)(3) Preparation and expiration dates. (c)(4) Other pertinent information required for proper use. This STANDARD is not met as evidenced by: Based on laboratory observation, review of the manufacturer's control package insert, and staff interview, the laboratory failed to label three of three control vials used to perform quality control on the Beckman Coulter DxH520 hematology analyzer with the revised expiration date on the survey date (08.19.2026). The findings include: 1. An observation of the laboratory on 08.19.2026 at 8:50 a.m. revealed a Beckman Coulter DxH520 hematology analyzer (Serial Number: BG060287) used for complete blood count (CBC) patient testing. Observation also revealed three levels of DxH500 Series hematology controls (lot numbers: abnormal low-352618711, normal-362618712, abnormal high-372618713) that lacked the revised expiration dates. 2. A review of the manufacturer's control package insert revealed that "16* Open Vial Days, *Assumes that the Instructions for Use section of the Consumable IFU/Setting Sheet is performed a maximum of 16 times within 16 days, provided they are handled properly." 3. An interview with the senior sight nursing supervisor on 08.19.2026 at 10:45 a.m. confirmed the above survey findings. Word key: IFU = Instructions for Use